BWT PERMAQ® compact 11-14

BWT PERMAQ® compact 15-18

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Table of contents

EN

Chapter 1: Introduction & safety	34
1.1 Abbreviations and subject index	34
1.2 Scope of delivery	35
1.3 Manufacturer	35
1.4 General comments	35
1.4.1 Reading of the manual (IOM)	35
1.4.2 Warranty provisions and disclaimer	36
1.4.3 Obligations of the operator	36
1.4.4 Licence conditions	36
1.4.5 Explanation of listed symbols	36
1.5. Operation and safety precaution advices	36
1.5.1 Intended use of device	36
1.5.2 Permitted mode of operation	36
1.5.3 Prohibited mode of operation	37
1.5.4 Qualified persons and users	37
1.5.5 Shutdown periods	37
1.6 Function of the RO unit	37
1.7 Preconditions for the installation	38
1.7.1 Positioning of the RO/installation requirements	38
1.7.2 Feed water requirements	38
1.7.3 Setting the correct operating pressure	38
Chapter 2: Assembly and installation	39
2.1 Preparation for assembly	39
Chapter 3: Operation and programming of the software	40
3.1 Quick operation guide	40
3.2 Function keys and display elements	40
3.3 Switching the device ON and OFF	41
3.4 Menu display of the activated system	41
3.5 System settings/functions	41
3.6 Main menu: SETTINGS (settings)	42
3.6.1 Menu: SETTINGS / GENERAL	42
3.6.2 Menu: SETTINGS / PRESSURE	44
3.6.3 Menu: SETTINGS / CONDUCTIVITY	46
3.6.4 Menu: SETTINGS / WCF SETUP	47
3.7 Main menu: INFO (information, statistics)	48
3.7.1 Menu: INFO / SYSTEM INFO	48
3.7.2 Menu: INFO / RESET	48
3.7.3 Menu: INFO / FACTORY DEFAULT SETTINGS	48
3.8 Main menu: INSTALLATION (installation by BWT service)	49
3.9 Main menu: MAINTENANCE (maintenance)	52
3.9.1 Menu: MAINTENANCE / MANUAL RINSING	52
3.9.2 Menu: MAINTENANCE / AUTOMATIC RINSING	52
3.10 Settings of the water quality using the blending valve	54

Table of contents

EN

Chapter 4: Fault finding	54
4.1 Trouble shooting guide.....	54
4.2 General check list for fault elimination.....	56
Chapter 5: Maintenance and service	57
5.1 Maintenance and wearing parts.....	57
5.2 Cleaning	58
5.3 Replacement of membrane element.....	58
5.4 Option: Cable kit for external alarm contact (OUTPUT)	58
5.5 Utilisation of the statistic timer	58
5.6 Disposal	58
Chapter 6: Technical data	59
6.1 Technical data BWT PERMAQ® compact 11-14	59
6.2 Technical data BWT PERMAQ® compact 15-18	60
Chapter 7: Documentation.....	61
7.1 Form for set values of BWT PERMAQ® compact 11-18	61
EC declaration of conformity.....	125

1.1 Abbreviations and subject index

Softening:

The water purification process removes the hardness from raw water. Hardness constituents are the portion of calcium and magnesium ions in the water.

Raw water:

Raw water (usually untreated drinking water) must often be pre-treated (usually softening) and is then useable for the desalination process in the RO device.

RO:

Abbreviation for Reverse Osmosis.

Permeate:

The largely desalinated "pure water" filtered in the RO membranes and generated by reverse osmosis. The characteristic value is the electric conductivity in $\mu\text{S}/\text{cm}$.

Concentrate:

Waste water led to outlet. This waste water contains the salts and minerals that have been removed from the raw water.

Membranes:

The "filter" of the device which is capable of desalinating the raw water by high pressure and flow.

SDI:

Abbreviation for "Silt Density Index". The "Silt Density Index" is an indicator of the organic contamination of raw water. The measurement system is a filtration process which determines the blocking tendency, measured during 15 min.

TDS:

Abbreviation for "Total Dissolved Solids" the total amount of dissolved salts, measured in mg/l .

Conductivity:

Electric conductivity value of the water, the smaller this measured value ($\mu\text{S}/\text{cm}$) the better the water quality.

IOM:

Abbreviation for "Installation and Operating Manual"

Permeate yield (WCF):

The ratio between the produced pure water (permeate) and the thus required amount of feed water (soft water) is expressed as permeate yield (WCF) or "Water Conversion Factor".

1.2 Scope of delivery

The reverse osmosis unit BWT PERMAQ® compact 11-18 consists of 3 core components:

- Measuring and control unit,
- High pressure pump,
- RO membranes and pressure vessels for online and offline operation.

Fig. 1: BWT PERMAQ® compact 11-14 - Front view of the device


Fig. 2: BWT PERMAQ® compact 15-18 - Rear view of the device

The scope of delivery of the RO device includes:

Microprocessor controller:

- 1 OLED display (monochrome)
- 2 Function button: ON/OFF and standby
- 3 Function button: Upward
- 4 Function button: Downward
- 5 Function button: Confirmation/OK
- 6 Function button: Return to previous menu
- 7 Duo-colour LED red/green

Connections and service components:

- 8 Permeate connection (Ø 8mm)
- 9 Concentrate connection (Ø 8mm)
- 10 Feed water connection (¾" or Ø 10mm)
- 11 Electrical wiring with grounded PE mains plug "CEE7"
- 12 Blending valve (adjustments only by BWT service)

Optionally available, Item no. 812850:

- 13 Installation kit for external alarm contact (OUTPUT), with M8 connector for the alarm output (DO).

1.3 Manufacturer

The RO devices BWT PERMAQ® compact 11 - 18 are manufactured in the name of:

BWT water + more GmbH

Walter-Simmer-Straße 4

AT-5310 Mondsee

Phone: +43/6232/5011-0

Fax: +43/6232/4058

E-Mail: info@bwt-filter.com

1.4 General comments

Following the instructions in this Installation and Operating Manual (IOM) helps the operator to run the reverse osmosis BWT PERMAQ® compact 11-18 reliably and economically. This Manual (IOM) is part of the device and must be constantly available at the place of operation for all assigned staff members.

1.4.1 Reading of the manual (IOM)

The staff must have read and understood this IOM prior to any work being carried out. A basic precondition for safe working is the adherence to all stated safety and operating instructions.

In addition, the local accident prevention provisions and the general safety provisions effective at the place of operation are applicable. The illustrations in these instructions serve the basic understanding and can deviate from the actual design of the device. Justified claims cannot be derived from the same.

1.4.2 Warranty provisions and disclaimer

Please observe: All of the information and instructions contained in this installation and operating manual were provided in respect of current standards, regulations, prior art and our long term experiences.

BWT excludes any responsibility for damages and for consequential loss due to:

- Non-compliance with the instructions in the manual
- Use not conforming to the intended purpose
- Improper or faulty installation
- Improper initial operation and operation, maintenance
- Use of non-permitted components as well as non-original parts
- Lack of examination in required service and maintenance work
- Damage due to unauthorized modification and technical manipulation

1.4.3 Obligations of the operator

- The installation and operating manual has always to be kept close to and accessible at the RO unit.
- The RO unit must only be used if it works in a technically proper manner and if the safety status is reliable.
- The instructions, notices and advices of the installation and operating manual have to be followed properly.

1.4.4 Licence conditions

The IOM is protected by copyright. Surrendering the manual to any third party, duplication of any kind and form – also in excerpts – as well as the utilisation and/or communication of the content are not permitted without the written consent of the manufacturer. Infringements obligate to pay compensation for damages. Further claims are reserved.

Note: Upon the acquisition of the reverse osmosis, the user receives an exclusive, non-transferrable right of utilisation in respect of the software installed by the manufacturer.

1.4.5 Explanation of listed symbols

This installation and operating manual uses **warnings/symbols!** The warnings/instructions are introduced by signal words to **underline risks**.

Please follow instructions and treat them with maximum importance to avoid any accidents and damages.

Danger: Caused by electric current or voltage!

Always consult a qualified electrician when working on places denoted by this symbol.

Attention: Dangerous spot!

Details or orders and prohibitions to avoid personal injury or extensive damage to property.

Please observe: Underlines useful recommendations and information for an efficient operation free of any interruptions.

Note: Additional information for the operator.

1.5 Operation and safety precaution advices

The following chapter provides an overview of all important operation and safety precaution advices to enable safe operation free of any interruptions. Despite all safety precaution measures, residual risks still remain, particularly in case of improper use and handling. A warranty claim will be denied, if the notices and advices of this installation and operating manual are not respected.

1.5.1 Intended use of device

The device is suitable for the desalination of water in drinking water quality up to **max. 30°C** and a raw water pressure of **max. 6.0bar** and can be operated either directly in front of the consumer in an **“online”** process or it can be connected **“offline”** to an atmospheric tank upstream of the user.

Observe: The feed water quality has to meet the requirements of the technical specifications (see Part 6.1 and Part 6.2) and must not be exceeded!

The device shall only be used for its intended purpose, in line with the installation and operating manual and in the environment to which it was intended and designed. Deviations are regarded as “non-intended purpose”.

The intended use also comprises the adherence to the operating, maintenance and upkeep conditions stipulated by the manufacturer.

Attention: The RO unit must be fed with cold water in drinking water quality.

- Every deviation from the intended purpose, e.g. deionisation of feed water of non-drinking water quality, can lead to irreversible damages or can cause undesired microbial contamination of the unit.

1.5.2 Permitted mode of operation

Observe: To protect the drinking water, observe the country-specific guidelines for drinking water installations in any work at the RO device.

- Prior to any maintenance work at the drinking water supply system, the device must be separated from the water supply system. Rinse the water pipe sufficiently before reconnecting the device.
- Interrupt the voltage supply of the device and of the terminal equipment prior to assembly (unplug mains plug).

Please observe! Improper and faulty installation of the RO unit can cause material damages.

- Please observe all applicable specific national guidelines and installation regulations (e.g. DIN 1988, EN 1717), general hygiene requirements and technical specifications in order to protect the drinking water.
- Unauthorized reconstruction work as well as technical modifications are prohibited.
- Please avoid any mechanical damages of the RO unit, otherwise the warranty will lapse.
- Please install a shutoff valve in front of the RO device.
- Please use **flexible tubes** according DVGW W543 requirements only.
- Chemicals, solvents and fumes must not come into contact with the RO unit.
- The installation site must be free of frost and exposure to sunlight.
- Do not use for example, **microbial contaminated feed water** or with **feed water of an uncertain quality or origin**.
- If permeate is used for **food preparation, clean and/or flush downstream users/machines before use**.
- **Avoid unnecessarily long storage times** of the equipment, in order to prevent the risk of microbial contamination.
- The demineralised water (permeate) **may not be used as drinking water**.

1.5.3 Prohibited mode of operation

Attention: DANGER due to use deviating from the intended purpose!

Any use beyond the intended purpose and/or any different use of the device can lead to dangerous situations.

Any claims because of damages resulting from any use not corresponding to the intended purpose are excluded.

Please observe:

NEVER operate the device with the housing covers removed.

The complete power supply unit must be replaced in the event of damage to the mains cable.

Please use:

- If needed, use protective clothing.
- **Do not use** any aggressive cleaning agents!

1.5.4 Qualified persons and users

Please observe:

Stipulate clear staff responsibilities for operation, set-up, maintenance, repair work!

For installation, commissioning and maintenance/service work, it is mandatory that it be effected by skilled and instructed person. The user has to be trained by a BWT organisation or other authorised person to operate the RO device.

- **Trained and instructed staff:** Was instructed about all possible dangers due to improper use.
- **Qualified persons:** Are able to install, commission and service an RO unit due to their qualification, knowledge and experience in current regulations.

1.5.5 Shutdown periods

The BWT PERMAQ® compact 11-18 unit is equipped with an programmable interval flush to prevent microbial contamination during long non-service periods. Nevertheless, it is recommended to implement the following measures in case of unfavourable cases.

- Remove the mains plug of the device and close the cut-off valve in the inlet of the reverse osmosis in case of longer breaks in operation.
- We recommend **discarding** the first **30 minutes of permeate** after each long idle periods of non use – e.g. weekends, holidays.
- Please carry out a replacement of the membrane elements, especially in unfavourable cases, e.g. longer down times in warm installation rooms.
- Please consider the individual manual of the external pre-filter.

1.6 Function of the RO unit

A semi-permeable RO membrane separates feed water of a high pressure (8...13 bar) into deionised (permeate) and into waste water (concentrate) with a high salinity.

The relation between the produced permeate and feed water is expressed as recovery or **WCF (%)**.

The RO unit can be operated either **(online)** (directly connected to the user) or **alternatively** connected with the universal level control **BWT levelmeter** (Order No. 812 854) to upfill the produced permeate into an atmospheric tank.

An integrated mixing valve enables blending of permeate and feed water in a precise manner to obtain a required product water conductivity depending on the application. Please take the recommended setting, see **Part 3.10 "Settings of the water quality with the blending valve"**.

The RO unit **starts** and **stops** according to the **values set for permeate pressure** in the control unit.

The control monitors several operating parameters, as soon as permeate is required, the pressure declines on the permeate side and the device starts automatically. As the turn-off pressure is reached, the device assumes the **"Ready"** stand-by mode.

If a failure occurs, the display shows the failure. If required, the device deactivates itself automatically. **Part 4.1** contains information concerning possible displays and their causes as well as **troubleshooting**.

1.7 Preconditions for the installation

1.7.1 Positioning of the RO/installation requirements

When installing the device, select a location where the device can easily be connected to the water supply network.

A connection to the sewage system and a separate mains socket (230V, 50Hz) must be nearby. Please connect the PE mains plug only to a grounded socket.

The rated **mains power** (see **technical data**) and the requisite **feed water pressure** must be present at all times.

If a pump station is needed for draining, please consider corrosion-resistant materials.

National guidelines and regulations:

Observe all applicable installation regulations, general guidelines, hygiene requirements, and technical specifications.

Frost protection and ambient temperature:

The installation site must be dry and free of frost and ensure the protection of the unit against chemicals, paint, solvents and fumes.

If the mains water is treated with oxidising disinfectants (chlorine, chlorine dioxide, etc.), an activated carbon filter is recommended to be placed in front of the unit. A **100µm particle filter** must always be provided by the customer.

Further pre-treatment is to be determined by the BWT service in dependence on the local feed water quality.

Quality of the pipeline supply:

Please observe: Only **corrosion-resistant material** must be used in the permeate area.

Electrical interference:

The emission of interference (voltage peaks, high-frequency electromagnetic fields, interference voltages, voltage fluctuations ...) by surrounding electrical systems may not exceed the maximum values specified in EN 61000-6-4.

1.7.2 Feed water requirements

The RO unit must only be fed with **cold drinking water (maximum 30°C)**, which meets the official requirements for potable water.

Data analysis of the feed water in your area:

Every deviation from the intended purpose, e.g. deionisation of non drinking water feed water can lead to irreversible damages to human and equipment (e.g. undesired microbial contamination of the RO unit).

Warning: Caused by unsuitable water quality.


- The limit values listed in Part 6.1 and 6.2 must not be exceeded.
- An exceeding of the limit-values by incompatible feed water might cause danger to unwanted sedimentation on the membranes. That influences negatively quality and performance.

1.7.3 Setting the correct operating pressure

A minimum operating pressure is required for the optimum functioning of the unit. Besides, the water pressure should not exceed the maximum permissible pressure.

Attention: The feed water inlet pressure must be between 2.0 and 6.0 bar!

If the pressure is more than 6.0bar, install a pressure reducing valve.

Please do not operate the unit if the pressure is below 2.0 bar.

Please avoid pressure variations of more than ± 0.5 bar.

- A ball valve is to be assembled on the inlet side of the device, so that the raw water supply may be interrupted for service purposes.
- Customers should install a feed water supply of at least 3/8" to avoid unreliable operation. An undersized inlet causes the danger of an interruption of the operation due to insufficient water pressure or insufficient flow, e.g. when rinsing the membranes.
- A pressure-reducing valve can decrease the flow rate negatively.

2.1 Preparation for assembly

Remove the device from its packaging and check the delivery to see that everything is there and that no damage was caused during transport.

Hydraulic installation:

Please respect the general installation instructions for the preparation of water installations as well as the general hygiene requirements.

- Please read/observe all applicable specifications, operating and safety instructions prior to the installation.
- Please use only approved flexible hoses, according to the DVGW W 543 requirements.
- Please adhere to all dimensions as well as bending radius in the assembly of the flexible tubes and connection sets.
- The BWT PERMAQ® compact 11-14 device line is to be installed and operated **horizontally** and the BWT PERMAQ® compact 15-18 device line **vertically**.
- The unit may not be connected to the water mains by rigid pipes.
- A hydrophore could be installed in the permeate line between the RO unit and the consumer, if depending on the application a short-term peak flow is required.

Connection to the water pipe:

- The device hoses must be installed in a flexible manner (without any tension) to maintain a safe operation.
- Check all water connections to ensure that they are tight and waterproof.
- Guide the flexible concentrate hose with **"free flow"** condition to the on-side connection of the drain pipe (with a slope of 1%) and fasten it. Note that the concentrate hose has been correctly connected. Please take care that the **"flexible" hose do not have any bends or mechanical reductions**.

Notes for initial operation:

- Connect the RO unit with the electrical power supply (230V, 50Hz). The wall socket has to be earthed.
 - Please consider the individual manual of the **external pre-filter**.
 - The **water hardness can vary** in different areas (soft water, hard water).
- ⇒ All hoses have been connected and are watertight.
- Please open the valve for the feed water supply.
 - Plug in the mains cable of the power supply.
- ⇒ The functions of the control unit and software configuration are described in **Part 3.1 - 3.9**.
- The setting of the permeate conductivity by blending is described in **Part 3.10**.
 - **Note: Discard the first permeate** produced during the first **30 minutes** after any new installation/first commissioning.

Note: Please run the RO unit a few days to reach the full performance (WCF and permeate quality).

Note: A drop of temperature by 1°C will reduce the permeate output of the membranes by approx. 3.0%.

Possible installation layout:

- 1 BWT PERMAQ® compact 11-18 reverse osmosis
- 2 External pre-filter (not included in the delivery)
- 3 Permeate outlet ready for the connection to the consumer/tank
- 4 Stop valve for feed water (not included in the delivery)
- 5 Mains plug, cable length 1.8m
- 6 Concentrate outlet
- 7 Optionally: External alarm contact (OUTPUT), alarm output


Fig. 3: Installation scheme / example

3.1 Quick operation guide

The BWT PERMAQ® compact 11-18 devices are controlled and parameterised using the control unit.

In this respect, the function keys listed below are available (see: [←], [▲], [▼], [✓], [ON/OFF]).

The illuminated alphanumeric display shows the main menu consisting of five symbols after the starting phase.

3.2 Function keys and display elements


Fig. 4: Display: Activated "Start" menu

[ON/OFF] key Switches the device from **standby** to ON/OFF and performs the reset of the device in case of alarms and/or failures.

[▲] (Arrow up) key Scrolls through the menu pages upward / right, increases the value of the parameter presently shown on the display.

[▼] (Arrow down) key Scrolls through the menu pages downward / left, decreases the value of the parameter presently shown on the display.

[✓] Confirmation key Confirmation of the parameter presently shown on the display.

[←] Return key Back to the previous menu or parameter-setting.

LED

LED for operating status and failure indication:

The status LED (green/red) shows the status of the device.

- **Green LED** for activated mains voltage (**LINE**)
- **Red LED** for failure messages (**FAILURE**)

Comments on the display:

Note: The display always switches OFF after **120s**.

Note: The menu **"FACTORY DEFAULT SETTINGS"** is secured by a password. Please contact the BWT service, in case of a change of the system settings.

Main menus:

Each of the **five main menus**, has its own symbols. The device is automatically switched to the menu item: **"Start"**.

- 1 Info:** *Information* (incl. timer and system reset)
- 2 Maintenance:** *Maintenance*
- 3 Start:** *Automatic mode* with target value input
- 4 Settings:** *Settings* (including general, pressure, conductivity, WCF target value input)
- 5 Installation:** *Installation menu* only to be use by BWT service


Fig. 5: Display: Main menu: "Ready: Press ✓ to start"

Note: To change into other menus, press the [▲] key or the [▼] key.

3.3 Switching the device ON and OFF

Switching the device ON:

The device is activated by pressing the **[ON/OFF]** key (approx. 2sec). Subsequently, the device is checked, rinsed and filled automatically.

Automatic start sequence of device:

BWT PERMAQ compact	Switching the RO unit ON
Autotest	Internal function check of sensors, electronics and pump
WATER FILL-UP	Filling the feed water into the RO unit
RINSING	Automatic rinsing of the RO unit

The RO unit shifts automatically to the **main menu "START"**.

• Main menu "Start":

Fig. 6: Display: Main menu: "Ready: Press ✓ to start"

Attention: Prior the initial start-up, ensure that the device has been checked for correct installation.

The system switches to the **"Starting"** status by tapping the **[✓]** **confirmation key** and is ready for operation.

Note: The production starts automatically as permeate is demanded (pressure drop on the permeate side).

- **Only if needed:** In the course of first commissioning the permeate conductivity should be set via the blending valve, see **Part 3.10**.

Switching the device "Standby":

Brief pressing of the **[ON/OFF]** key, the device switches to standby. The device is immediately available for operation after reactivation.

Switching the device OFF:

Deactivate the device completely by pressing the **[ON/OFF]** key for **about 2s** on the control unit.

3.4 Menu display of the activated system

PERMEATE	1975 ml/min	COND	12µS
CONCENTRATE	1190 ml/min	TEMP	17°C
WCF SET POINT	62 %	WCF	62%
PUMP PRESS	13.4 bar	= Working	
PERM PRESS	1.1 bar		

Starting
Working
Ready

Fig. 7: Display: Activated menu "Start" at process step: **Working**

Status after activation the "Start" menu:

Starting:	RO unit starts
Working:	Permeate production is activated
Ready:	Ready for a new consumption demand

Note: The display always switches OFF after about **120s**. A brief pressing of any key re-activates the display.

3.5 System settings/functions

Note: All of the settings described in **Part 3.8** can only be accessed and activated in the **"Ready"** mode and the **"START"** menu of the device.

During operation mode **"Working"** are any parameter settings NOT possible.

Visualisation of the menu screens/displays:

Note: The differing representation of the software diagrams exclusively serves a better overview and is subsequently described.

Bold black writing is subsequently used to highlight respectively **activated menu items or data input**. The **inactive menu** items are represented in **dark grey**.

• GENERAL	ACTIVE
• PRESSURE	INACTIVE

Blue fields = Menu item of a main menu, that contains further sub-menus or parameter inputs.

White fields, in grey writing = inactive menu.

Black writing, bold = active menu line or input value.

Dark grey writing = inactive menu line or input value.

3.6 Main menu: SETTINGS

General settings:

⇒ The main selection menu appears with "Start" being activated.


Fig. 8: Main selection menu: "Start", "Settings"

The "Settings" selection menu appears:

⇒ The menu item "● GENERAL" is active.

● GENERAL	General device settings (ACTIVE)
● PRESSURE	Pressure of the permeate
● CONDUCTIVITY	Conductivity value (µS/cm)
● WCF SETUP	Permeate output/yield WCF menu

Fig. 9: Display: Activated "Settings" menu

The first menu level of the "Settings" menu consists of four selectable menu items in which the subsequent settings may be effected:

GENERAL:

- General device settings (date, time, language)

PRESSURE:

- Set point and limit values for pressure

CONDUCTIVITY:

- Set point and limit values for conductivity

WCF SETUP:

- Set point value for recovery (WCF), efficiency in %

3.6.1 Menu: SETTINGS / GENERAL

Perform the general settings of the selection menu "Settings / GENERAL" in accordance with Fig. 10.

▶ DATE	Date (ACTIVE)
▶ TIME	Time
▶ LANGUAGE	Language

Quick operation overview, applies to all menus:

[✓] key ... for activation of all menus.

[←] key ... to return into the main menu.

Note: A brief pressing of any key re-activates the display.

Note: At present, the software is only available in the English language.


Fig. 10: Programming procedure: "GENERAL"

3.6.2 Menu: SETTINGS / PRESSURE

Settings for permeate pressure parameters:

⇒ The main selection menu appears with **"Start"** being activated.


Fig. 11: Main selection menu: "Start", "Settings"

1. Menu level (selection menu "Settings")

● GENERAL	General device settings
● PRESSURE	Pressure of the permeate (ACTIVE)
● CONDUCTIVITY	Conductivity value ($\mu\text{S}/\text{cm}$)
● WCF SETUP	Permeate output/yield WCF menu

Fig. 12: Menu Level 1: "Settings"

Perform the general settings of the **"Settings / PRESSURE"** selection menu in accordance with **Fig. 14**.

The recommended target values of all settings are documented in the form "Form for set values of BWT PERMAQ® compact 11-18" (see Part 7.1).

2. Menu level (adjustable parameters)

▶ START PRESSURE	Start pressure (ACTIVE)
▶ STOP PRESSURE	Stop pressure
▶ STOP TIME DELAY	Time delay of the pump

Fig. 13: Menu Level 2: "Settings/PRESSURE"


Fig. 14: Programming procedure: "PRESSURE" (permeate pressure)

3.6.3 Menu: SETTINGS / CONDUCTIVITY

Settings for permeate conductivity:

⇒ The main selection menu appears with "Start" being activated.


Fig. 15: Main selection menu: "Start", "Settings"

⇒ The menu item "● CONDUCTIVITY" is active.

1. Menu level (selection menu "SETTINGS")

● GENERAL	General device settings
● PRESSURE	Pressure of the permeate
● CONDUCTIVITY	Conductivity value (ACTIVE)
● WCF SETUP	Permeate output, WCF menu

Fig. 16: Menu Level 1: "Settings"

Press the [✓] key again to set the limits for the alarm and the activation/deactivation of the conductivity alarm.

2. Menu level (adjustable parameters)

▶ CONDUCTIVITY LIMIT	Conductivity limit value (ACTIVE)
▶ CONDUCTIVITY ALARM	Conductivity alarm

Fig. 17: Menu Level 2: "Settings/CONDUCTIVITY"

Note: In case you activate the "STOP" setting, you have to deactivate and then activate the device again.

Note: It is not possible to set the last digit "100" at first. Please observe the right input order.

Perform the settings of the "CONDUCTIVITY" selection menu in accordance with Fig. 18.


Fig. 18: Programming procedure: "CONDUCTIVITY"

3.6.4 Menu: SETTINGS / WCF SETUP

• Settings for recovery (WCF):

⇒ The main selection menu appears with "Start" being activated.


Fig. 19: Main selection menu: "Start", "Settings"

⇒ The menu item "● WCF SETUP" is active.

1. Menu level (selection menu "SETTINGS")

● GENERAL	General device settings
● PRESSURE	Pressure of the permeate
● CONDUCTIVITY	Conductivity value ($\mu\text{S}/\text{cm}$)
● WCF SETUP	Permeate output WCF (ACTIVE)

Fig. 20: Menu Level 1: "Settings"

Press again the [✓] key to adjust the set points of the WCF menu.

Please observe: It is only possible to enter values within the programmed range of values.

Perform the settings of the "WCF SETUP" selection menu in accordance with Fig. 21.


Fig. 21: Programming procedure: "WCF SETUP"

3.7 Menu: INFO (information, statistics)

⇒ The main selection menu appears with "Start" being activated.


Fig. 22: Main selection menu: "Start", "Info"

⇒ The menu item "► SYSTEM INFO" is active.

SYSTEM INFO:

- System information (date, time, language)

3.7.1 Menu: INFO / SYSTEM INFO

In the "Info" selection menu, you can both read system information/device data and counters and reload factory parameters.

The recommended target values of all settings are documented in the form "Form for set values of BWT PERMAQ® compact 11-18" (see part 7.1).

1. Menu level (selection menu "Info")

► SYSTEM INFO	System information (ACTIVE)
► RESET	Reset of statistic timer (BWT only)
► FACTORY DEFAULT SETTINGS	Save factory settings as standard (BWT only)

Fig. 23: Menu Level 1: "Info"

2. Menu Level (password-protected parameters)

► PASSWORD	System password
*****	Enter the password:

PASSWORD ENTRY: (down arrow) (down arrow) (up arrow) (up arrow) (up arrow) (down arrow)

BWT Service: Please select and set each digit individually. Firstly, by selection with the [▲] (arrow up) and [▼] (arrow down) keys and, secondly, by confirming each digit.


Fig. 24: Programming procedure: "INFO" (system information)

3.7.2 Menu: INFO / RESET

⇒ The second menu item "► RESET" is active.

- BWT service **reset** of the **statistic timer** (deletes all values).

► PASSWORD	System password
*****	Enter the password:

3.7.3 Menu: INFO / FACTORY DEFAULT SETTINGS

⇒ The 3rd menu item "► FACTORY DEFAULT SETTINGS" is active.

- Save parameters as new standard factory settings.

► PASSWORD	System password
*****	Enter the password:

3.8 Menu: INSTALLATION

⇒ The main selection menu appears with "Start" being activated.


Fig. 25: Main selection menu: "Start", "Installation"

⇒ The menu item "► WCF MIN" is active.

The following operational parameters may be set in the "INSTALLATION" selection menu (Fig. 30):

► WCF MIN	Minimum value for permeate yield
► WCF MAX	Maximum value for permeate yield
► PRESSURE STOP MAX	Deactivation pressure RO
► PRESSURE START MAX	Activation pressure RO
► STOP TIME DELAY	Stop delay time (HP-pump)
► MAX PUMP PRESSURE	Maximum pump pressure
► RINSE TIME	Step time: Rinsing
► WATER FILL UP TIME	Step time: For feed water filling

Fig. 26: Menu Level 1: "INSTALLATION"

You effect settings in the "INSTALLATION" menu according to the following overview menu structure of Fig. 27, Fig. 28, Fig. 29.

Note: The settings in the "INSTALLATION" menu may only be effected by qualified and instructed persons.

Please observe: The setup menu is **not password-protected**. Please **enter** the **newly set parameters** in the form of **Part 7.1**.

The adjustment of "WCF MIN" and "WCF MIN" is effected according to BWT specifications in accordance with the available water quality.


Fig. 27: Programming procedure: INSTALLATION 1/3


Fig. 28: Programming procedure: INSTALLATION 2/3


Fig. 29: Programming procedure: INSTALLATION 3/3

3.9 Main menu: MAINTENANCE

⇒ The main selection menu appears with **"Start"** being activated.


Fig. 30: Main selection menu: "Start", "Maintenance"

3.9.1 Menu: MAINTENANCE / MANUAL RINSING

⇒ The first menu item **"► MANUAL RINSE"** is active.

► MANUAL RINSING	<i>Mode for manual flush (ACTIVE)</i>
► AUTOMATIC RINSING	<i>Mode for automatic flush</i>
► CLEANING	<i>Cleaning deactivated</i>
► DISINFECTION	<i>Disinfection deactivated</i>

Fig. 31: Menu Level 1: "Maintenance"


Adjustment for "MANUAL RINSING":
ATTENTION: The manual flushing begins immediately.

⇒ The RO unit rinses manually.

3.9.2 Menu: MAINTENANCE / AUTOMATIC RINSING

In this menu, **"AUTOMATIC RINSING"**/interval rinsing at a set time is programmed.

Please perform the setting in the menu **"AUTOMATIC RINSING"** according to the **Fig. 32**.

Document the parameters after you have set all of the data at the device.

Please use the template at the end of the operating instructions.


Fig. 32: Programming procedure: "MAINTENANCE"

3.9 Main menu: MAINTENANCE

⇒ The main selection menu appears with **"Start"** being activated.


Fig. 30: Main selection menu: "Start", "Maintenance"

3.9.1 Menu: MAINTENANCE / MANUAL RINSING

⇒ The first menu item **"► MANUAL RINSE"** is active.

► MANUAL RINSING	<i>Mode for manual flush (ACTIVE)</i>
► AUTOMATIC RINSING	<i>Mode for automatic flush</i>
► CLEANING	<i>Cleaning deactivated</i>
► DISINFECTION	<i>Disinfection deactivated</i>

Fig. 31: Menu Level 1: "Maintenance"

Adjustment for "MANUAL RINSING":
ATTENTION: The manual flushing begins immediately.

⇒ The RO unit rinses manually.

3.9.2 Menu: MAINTENANCE / AUTOMATIC RINSING

In this menu, **"AUTOMATIC RINSING"**/interval rinsing at a set time is programmed.

Please perform the setting in the menu **"AUTOMATIC RINSING"** according to the **Fig. 32**.

Document the parameters after you have set all of the data at the device.

Please use the template at the end of the operating instructions.


Fig. 32: Programming procedure: "MAINTENANCE"


Fig. 33: Blending valve: "H₂O-Mix" - Rear view of the device

3.10 Settings of the water quality using the blending valve

The permeate conductivity is increased by the addition of feed water via the **integrated blending valve (12)**.

This blending is required depending on the intended use.

The screw is on rear side of the device, **see Fig. 33**.

The setting is changed using a **hexagon socket wrench (3mm)**.

The conductivity of the permeate can be adjusted:

Rotation of the control valve to the left (+)

increases the value of the conductivity.

Rotation of the control valve to the right (-)

decreases the value of the conductivity.

The current conductivity value automatically appears on the display in $\mu\text{S}/\text{cm}$ during water intake.

4.1 Trouble shooting guide

Please observe: If your reverse osmosis is not performing as it should, please run through the **alarm messages** below and the **check list (Part 4.2)** on **page 26**. In the event of a fault, repair work may only be carried out by an expert (SERVICE technician).

Alarm / malfunction list:	Possible reasons:	Acknowledgement:
FAULT 1 BUS COMMUNICATION ERROR	— <i>FAULT 1 ... in operation</i>	
FAULT 1 BUS COMMUNICATION ERROR 	— <i>FAULT 1 ... during starting</i>	
1: "BUS COMMUNICATION ERROR" • Bus communication error	• Cable connection between measuring board and control interrupted or defective	⇒ Call SERVICE team ⇒ Repair defect parts
FAULT 2 PUMP FUSE	— <i>FAULT 2 ... in operation</i>	
FAULT 2 PUMP FUSE 	— <i>FAULT 2 ... during starting</i>	
2: "PUMP FUSE" • Fuse error of the high-pressure pump	• Motor fault, check operating parameters • Pump is mechanically blocked or damaged	⇒ Call SERVICE team ⇒ Replacement of pump motor

4.1 Trouble shooting guide

 **Please observe:** If your reverse osmosis is not performing as it should, please run through the **alarm messages** below and the **check list (Part 4.2)** on the **next page**. In the event of a fault, repair work may only be carried out by an expert (SERVICE technician).

Alarm / malfunction list:	Possible reasons:	Acknowledgement:
FAULT 3 PUMP PRESSURE TRANSMITTER FAILURE FAULT 3 PUMP PRESSURE TRANSMITTER FAILURE 	— <i>FAULT 3 ... in operation</i> — <i>FAULT 3 ... during starting</i>	
3: "PUMP PRESSURE TRANSMITTER FAILURE" • Failure of the pump pressure transmitter	• Defective pressure transmitter or loose contacts	⇒ Call SERVICE team ⇒ Replacement of pressure transmitter
FAULT 4 MAX PRESSURE TRANSMITTER FAILURE FAULT 4 MAX PRESSURE TRANSMITTER FAILURE 	— <i>FAULT 4 ... in operation</i> — <i>FAULT 4 ... during starting</i>	
4: "MAX PRESSURE TRANSMITTER FAILURE" • Failure of the maximum pressure of the pressure sensor	• Pressure transmitter is defective • Faulty wiring of the pressure transmitter	⇒ Call SERVICE team ⇒ Replacement of pressure transmitter
PUMP PRESSURE TO HIGH Reduce it and press V to restart 5: "PUMP PRESSURE TOO HIGH" • The pump pressure is too high TEMPERATURE PROBE FAILURE 	— <i>FAULT 5 ... only in operation</i> • Effective pump pressure is larger than the limit set point • Pressure switch defective — <i>FAULT 6 ... only during starting</i>	⇒ Call SERVICE team ⇒ Open the bypass at the pump head to reduce the pressure and adapt the set point (max. 13bar) ⇒ Replacement of pressure transmitter
6: "TEMPERATURE SENSOR FAILURE" • Failure of the temperature sensor	• Temperature sensor defective	⇒ Call SERVICE team ⇒ Replacement of temperature sensor

4.1 Trouble shooting guide

 **Please observe:** If your reverse osmosis is not performing as it should, please run through the **alarm messages** below and the **check list (Part 4.2)**. In the event of a fault, repair work may only be carried out by an expert (SERVICE technician).

Alarm / malfunction list:	Possible reasons:	Acknowledgement:
 WATER LEAKAGE Disconnect power supply Close water feed	<i>FAULT 7 ... only in operation</i>	
7: "WATER LEAKAGE" • System water leakage	• Water leaks from the device	⇒ Call SERVICE team ⇒ Separate the unit from the power-supply and the water network
 CONDUCTIVITY > xx uS Press V to start	<i>FAULT 8 ... only in operation</i>	
8: "CONDUCTIVITY > xx µS" • Conductivity too high, the limit value has been exceeded (µS/cm)	• To technical data unconform feed water • Damaged conductivity sensor	⇒ Call SERVICE team ⇒ Please check the pre-treatment of the feed water.

4.2 General check list for fault elimination

Check list:	Possible reasons:	Rectification:
Problem: No water at the inlet connection?	• Shut-off valve of supply tap is closed, • Interruption of water supply	⇒ Open the main valve or water supply.
Problem: Permeate output too low? Problem: Is the line pressure within the specification of the device?	• Water inlet pressure too high • Water inlet pressure too low	⇒ Check that the water supply pressure is between 2.0 to 6.0 bar.
Problem: Only in case of softwater inlet: The feed water still remains hard.	• External water softener supplies hard water	⇒ Check the water softener and fill the brine cabinet of the water softener with salt.
Problem: Device flushes at the wrong time.	• Is the present time correct?	⇒ Reset the present time.

Note: If the problem persists, please call your local dealer or after sales service team.

5.1 Maintenance and wearing parts

You have purchased a product that is durable and easy to service.

However, all technical equipment requires regular servicing in order to guarantee optimal functionality.

If the product malfunctions during the warranty period, contact your contract partner, the installation company, and quote the unit type and production number (see technical specifications or the type plate on the unit).

Please exchange the wearing parts in the specified maintenance intervals:

Please observe!

- In order to safeguard unobjectionable functioning and optimum water quality, **maintenance** is to be performed at **regular intervals** by our service technicians, **at least however, once a year**.

- ✓ **Consider as well the instruction manual of the external installed pre-filter.**
- ✓ Prior to working on electric components and a opened housing, **ensure carefully that the device has been deactivated by unplugging the mains plug and therefore has zero potential.**
- ✓ **Any maintenance work should include a check of the connecting lines and the housing.**

Note: The replacement of the wearing parts must be effected by **qualified compartment staff**. (Plumber or customer BWT service).

Exchange of wearing parts:

Maintenance work:

- ✓ General visual inspection
- ✓ Check of tightness
- ✓ Cleaning with a moist cloth
- ✓ Conductivity (value displayed by RO)
- Conductivity (value of measuring device)
- Exchange of the external pre-filter insert (100µm)
- Exchange RO membrane elements
- Exchange of high-pressure pump head
- Exchange of connections/pressure hoses
- Water hardness test

Responsible:

- Customer
- Customer
- Customer
- Customer
- Customer, service
- Customer, service
- Service
- Service
- Service
- Service

Recommended maintenance interval:

- Weekly
- Weekly
- As required
- Weekly
- At least, once a year
- At least, 2 times a year
- At least, once a year
- Every 2 years
- Every 5 years
- At least, once a year

Please observe: According to BGV A2 (VBG4), the electrical safety check is to be conducted every 4 years.

The BWT PERMAQ® compact 11-18 devices have been designed in compliance with the "Pressure Equipment Directive" 97/23/EC of 29th May 2007. Therefore, it meets the requirements of Article 3, Section 3, and is consequently manufactured and assembled in accordance with applicable good engineering practice.

The BWT PERMAQ® compact 11-18 devices do not receive a CE mark in accordance with Article 15 of Directive 97/23/ EC. However, the attached EC declaration of incorporation is applicable.

Note: To ensure reliable production and to meet the requirements of our customers or to extend the warranty, our company offers the conclusion of a **BWT service contract**.

- ✓ **BWT service or a contract guarantees the professional, regular examination** in which all important components are cleaned, checked and in which wearing parts are replaced.
- ✓ Please ask about our services.

5.2 Cleaning

Clean the reverse osmosis device only with a ✓ **moist cloth** and **mild detergents**.

✓ **Use cloths without any fibres!**

✓ **Please do not use:**

- Bleaches,
- Solvents,
- Alcohol,

to protect the surface of the device.

5.3 Replacement of membrane element

In dependence of the feed water quality and the pre-treatment salts (predominantly calcium) can be deposited on the membrane and impair the performance (volume flow) of the membrane and the quality of the permeate.

If the permeate volume flow declines or the conductivity in the permeate rises, the membrane element has to be exchanged. **However, this is required latest after 12 months in any case.**

➤ **The exchange should be executed by the BWT service.**

Please record in case of a membrane exchange:

- ✓ **1.** Date of replacement of the new membrane element:
- ✓ **2.** New output of the RO unit (display):
- ✓ **3.** Water quality (µS/cm):
- ✓ **4.** Operating pressure of the device:
- ✓ **5.** Feed water temperature:
- ✓ **6.** Feed water pressure:
- ✓ **7.** Measurement of water hardness °dH (with hardness test set):

5.4 Option: Cable kit for external alarm contact (OUTPUT)

The optional **“Cable kit for external alarm contact (OUTPUT)”** enables the use of the potential-free contact for the alarm output. All RO alarms can be transferred as an common alarm to an external control.

The **“Cable kit for external alarm contact (OUTPUT)”** is either already preinstalled or upgradable as an additional install option.

Cable kit for external alarm contact (OUTPUT),
Order no.: 823428

5.5 Utilisation of the statistic timer

To activate the statistics timer tap the menu:
“INFO” / “SYSTEM INFO”.

permeate 03003 lit	avg wcf xx %
pump worked for	x hours
pump on	x times
Software ver 2.2	
Service maintenance 364 days	

The following maintenance informations are accessible:

permeate *** lit:**

- Sum of permeate in liters,

avg wcf xx %:

- Average recovery WCF in %,

pump worked for x times:

- Operating hours of the pump,

pump on x times:

- Counter of all pump starts,

Software ver 2.70:

- Software version,

Service maintenance * days:**

- Readout of the next Service, if the counter shows “0 days”.

5.6 Disposal

Procedure:

The BWT PERMAQ® compact 11-18 devices consists of various materials which need to be disposed of properly.

- ✓ Please order the manufacturer customer service for an expert and environmentally compliant disposal. Please do not throw used batteries into household waste.

- ✓ Disposal of any electrical parts should only be carried out at authorised WEEE recycling centres (EN 2002/96/EC).

6.1 Technical data BWT PERMAQ® compact 11 - 14


Technical data BWT PERMAQ® compact 11 - 14				
BWT PERMAQ® compact (horizontal)		11	12	14
Nominal capacity ^{*1)} (output flow rate)	l/h	50	75	100
Salt retention rate	%	> 95	> 95	> 95
Permeate output/yield (WCF) ^{*2)}	%	30 ... 75	30 ... 75	30 ... 75
Feed water flow (inlet)	l/h	≥ 67 ... 167	≥ 100 ... 250	≥ 167 ... 333
Concentrate flow (drain)	l/h	≥ 17 ... 117	≥ 25 ... 175	≥ 67 ... 233
Feed water pressure	bar	2.0 ... 6.0		
Feed water, ambient temperature (min./max.)	°C	5 ... 30 / 5 ... 40		
Iron and manganese (Fe+Mn)	mg/l	< 0.05		
Silica (SiO ₂)	mg/l	< 15		
Salt level, Total Dissolved Solids (TDS)	mg/l	< 500		
Silt Density Index (SDI)	%/min	< 3		
Oxidants	mg/l	< 0.05		
Protection class	IP	55	55	55
Electrical connection / fuse protection	V/Hz/A	230 / 50 / 10		
Electric power consumption (operation / standby)	W	460 / < 3	460 / < 3	460 / < 3
Type of valve		230VAC		
Standard of connector (grounded PE mains plug)		Three-pin "CEE7"		
Feed water, permeate, concentrate connection	inch/mm/mm	3/4" or 10mm / 8 / 8		
Dimensions: Width, depth, height (W x D x H)	mm	390 x 400 x 150	390 x 540 x 150	390 x 540 x 150
Weight	kg	18	20	20
Production and Item No. (BWT AT)		812811	812812	812814

Note: The nominal flow is calculated with following drinking water conditions:

^{*1)} Water inlet temperature 15°C, salinity as NaCl max. 500 mg/l, Primary pressure inlet 3.5bar and a permeate yield of 30%. During operation, the actual nominal flow may deviate slightly from the flow indicated in the table due to fluctuations of the inlet water quality, the flow pressure and the flow continuity.

^{*2)} For normal operation the manufacturer recommends the inlet of softened feed water with a hardness lower than 14°dH.

6.2 Technical data BWT PERMAQ® compact 15 - 18


Technical data BWT PERMAQ® compact 15 - 18				
BWT PERMAQ® compact (vertical)		15	16	18
Nominal capacity ^{*1)} (output flow rate)	l/h	50	75	100
Salt retention rate	%	> 95	> 95	> 95
Permeate output/yield (WCF) ^{*2)}	%	30 ... 75	30 ... 75	30 ... 75
Feed water flow (inlet)	l/h	≥ 67 ... 167	≥ 100 ... 250	≥ 167 ... 333
Concentrate flow (drain)	l/h	≥ 17 ... 117	≥ 25 ... 175	≥ 67 ... 233
Feed water pressure	bar	2.0 ... 6.0		
Feed water, ambient temperature (min./max.)	°C	5 ... 30 / 5 ... 40		
Iron and manganese (Fe+Mn)	mg/l	< 0.05		
Silica (SiO ₂)	mg/l	< 15		
Salt level, Total Dissolved Solids (TDS)	mg/l	< 500		
Silt Density Index (SDI)	%/min	< 3		
Oxidants	mg/l	< 0.05		
Protection class	IP	55	55	55
Electrical connection / fuse protection	V/Hz/A	230/50/10		
Electric power consumption (operation / standby)	W	460 / < 3	460 / < 3	460 / < 3
Type of valve		230VAC		
Standard of connector (grounded PE mains plug)		Three-pin "CEE7"		
Feed water, permeate, concentrate connection	inch/mm/mm	3/4" or 10mm / 8 / 8		
Dimensions: Width, depth, height (W x D x H)	mm	150 x 400 x 390	150 x 540 x 390	150 x 540 x 390
Weight	kg	18	20	20
Production and Item No. (BWT AT)		812815	812816	812818

Note: The nominal flow is calculated with following drinking water conditions:

^{*1)} Water inlet temperature 15°C, salinity as NaCl max. 500 mg/l, Primary pressure inlet 3.5bar and a permeate yield of 30%. During operation, the actual nominal flow may deviate slightly from the flow indicated in the table due to fluctuations of the inlet water quality, the flow pressure and the flow continuity.

^{*2)} For normal operation the manufacturer recommends the inlet of softened feed water with a hardness lower than 14°dH.

7.1 Form for set values of BWT PERMAQ® compact 11-18

Installation date:	 / /	<u>Fitter:</u>
Customer:		
Model:		
OBSERVE: THE SETTINGS 1-8 MAY BE CHANGED BY THE CUSTOMER & THE SETTINGS 9-19 ONLY BY THE BWT SERVICE.		
Parameter:	Value set at installation:	Description:
DATE:	[Day / month / year]	Current date
TIME:	[Hours : minutes]	Current time
LANGUAGE:	Language selection	Active language: English
START PRESSURE:	bar	Pmin = 2.0 bar
STOP PRESSURE:	bar	Pmax = 3.0 bar
STOP TIME DELAY:	seconds	Default setting: 10 s
CONDUCTIVITY LIMIT:	µS/cm	Default setting: 500 µS/cm
CONDUCTIVITY ALARM:	DISABLED / <u>ALARM</u>	Default setting: ALARM enabled
<u>WCF set point:</u>	Range: 20 ... 79 %	Default setting: 50 %
<u>WCF MIN:</u>	Range: 20 ... 59 %	Default setting: 40 %
<u>WCF MAX:</u>	Range: 59 ... 79 %	Default setting: 60 %
<u>PRESSURE START MAX:</u>	Range: 1.0 ... 7.0 bar	Default setting: 2.0 bar
<u>PRESSURE STOP MAX:</u>	Range: 1.0 ... 7.0 bar	Default setting: 4.0 bar
<u>STOP TIME DELAY:</u>	Range: 0 ... 19 s	Default setting: 10 s
<u>MAX PUMP PRESSURE:</u>	Range: 8 ... 13 bar	Default setting: 9 bar
<u>RINSING TIME:</u>	Range: 0 ... 599 s	Default setting: 10 s
<u>WATER FILLING TIME:</u>	Range: 0 ... 99 s	Default setting: 10 s
<u>AUTOMATIC RINSING:</u>	ENABLE / DISABLE	Default setting: ENABLE
<u>RINSING HOUR:</u>	TIME: 23:30 - 3:00	Default setting: 02:00 a.m.

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CE Declaration of conformity

EG Konformitätserklärung / Certificat de conformité CE

The company **BWT water + more GmbH** declares, that **reverse osmosis device** with the following specifications:

Trade name of product:	Model:	Order No.
• BWT PERMAQ® compact	11	812811
• BWT PERMAQ® compact	12	812812
• BWT PERMAQ® compact	14	812814
• BWT PERMAQ® compact	15	812815
• BWT PERMAQ® compact	16	812816
• BWT PERMAQ® compact	18	812818

with a **serial number higher than:** see rating plate & technical specifications

and with a order No.: see rating plate & technical specifications

have been **designed, manufactured** and **assembled** according to the following **EC Directives (guidelines):**

2004/108/EC Guideline for electromagnetic compatibility (EMC)

2006/95/EC "Low Voltage Directive" 2006/95/EC

the following harmonised guidelines were applied:

EN 61000-6-1, EN 61000-6-3, EN 60335

Manufacturer: BWT water + more GmbH, Walter-Simmer-Straße 4,
A-5310 Mondsee, Phone: + 43 (0) 6232 5011-0

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[illegible]

BWT – The Company

The Best Water Technology Group is Europe's leading water technology company with 2,800 employees and a vast network of partner companies, service staff, installers, planners, architects and hygiene experts. Our R&D teams apply state-of-the-art methods and use the latest processes and materials to create products that are both economical and green. One of the key objectives is to reduce the resource and energy consumption levels of our products, resulting in a reduction in CO₂ emissions.

Practically everywhere where water is involved, BWT's pioneering products have proved their worth millions of times over; where domestic water enters a building ('point of entry') and at its tapping point ('point of use'), for the treatment of seawater, drinking water, mineral water, ultrapure water for pharmaceuticals, water for swimming pools, heating and process water, boiler water, cooling water and water for air-conditioning systems. Our numerous innovations enable our customers to enjoy the highest levels of safety, hygiene and health in their daily use

of water – the precious elixir of life. Among these are: SEP-TRON®, the world's first Electrodeionisation module (EDI) with spiral winding, the MDA (manganese oxide activation) method for effective removal of manganese, AQA total bipolar technology for chemical-free limescale protection, SANISAL – the world's first regeneration salt for softening systems that also disinfects at the same time, and the revolutionary new Mg²⁺ technology for better tasting filtered water, coffee and tea. BWT, with its unique, high-performance membranes for fuel cells and batteries, is setting the 21st century standard for the supply of clean energy.

BWT – For You and Planet Blue signifies our mission to take ecological, economic and social responsibility, and to provide our customers with the best products, systems, technologies and services in all areas of water treatment, while making a valuable contribution to the preservation of our blue planet's global resources.

BWT water + more GmbH

Walter-Simmer-Straße 4
A-5310 Mondsee
Phone: +43 / 6232 / 5011 0
Fax: +43 / 6232 / 4058
E-Mail: office@bwt.at

Händlerstempel (DE) / Stamp of dealer (EN)
Cachet du revendeur (FR) / Timbro del rivenditore (IT)

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Indirizzo del produttore: BWT water + more GmbH, Walter-Simmer-Straße 4,